



GREER™
Kimball® Health

Greer is a modern and modest, full function patient recliner built with quality and value in mind. The small-scale design works perfectly for today's small footprint spaces while still providing comfort for all body shapes and sizes.

Infinite Back Release

- Pull lever out and use body weight to push back to desired recline, then release lever
- To return to upright position, pull lever out and lean forward or assist patient to upright position

Footrest Release Lever

- Pull back to lift footrest
- Push forward to return to original position

Trendelenburg Paddle

- Use footrest release lever to bring footrest completely out
- Use infinite back release to bring patient to lay-flat position
- Engage trendelenburg paddle with foot and guide patient to trendelenburg position
- Re-engage paddle and guide patient back to standard reclined position




Central Locking Caster

- Push red paddle down to lock all 4 casters
- Push green paddle to neutral position to engage non-directional movement
- Push green paddle forward to engage directional movement



Fold Down Tray

- Manually lift tray to upright position
 - Compress levers below table and guide table back down to closed position
- *Not intended to be used for patient transport**



Pull-Out Footrest

- Manually pull out for additional patient support
- *Not intended to be used for patient transport**



Drop Arm

- To release, push lever in at back of arm and guide arm to dropped position
- To close, lift arm to starting position, the lever will re-engage when put back to starting position